

著書

- Tanimizu, N. “Development and Regeneration of Intrahepatic Bile Ducts” Leon V. Berhardt, Ed. In *Advances in Medicine and Biology*. Volume 130, Nova Science Publishers, Chapter 6, pp. 187-216 (2018)
- Tanimizu N and Mitaka T. “Plasticity of Liver Epithelial Cells in Healthy and Injured Livers. In *Stem Cells and Cancer in Hepatology 1st Edition: From the Essentials to Application*. Ed. Yun-wen Zheng. Academic Press, p35-54 (382 pages) (2018) doi.org/10.1016/B978-0-12-812301-0.00003-7
- Mizuguchi T, Mitaka T, Hirata K. Role of Branched Chain Amino Acids in Cellular and Organ Damage: The prognostic significance of the preoperative branched chain amino acid to tyrosine ratio. In Rajendram R, Preedy VR, Patel VB, Eds. In *Branched Chain Amino Acids in Clinical Nutrition, Vol II*. London: Humana Press; 2015, pp65-77 (総ページ数 320)

総説

1. Tanimizu N, Mitaka T. Which is better source for functional hepatocytes? (Editorial) *Stem Cell Invest*, Feb 16; 4:12 (2017) doi: 10.21037/sci.2017.02.08
- 2.
3. Tanimizu N, Mitaka T. Epithelial morphogenesis during liver development. *Cold Spring Harb Perspect Biol*, 9(8). pii: a027862 (2017) doi: 10.1101/cshperspect.a027862
- 4.
5. Tanimizu N, Mitaka T. Morphogenesis of Liver Epithelial Cells. *Hepatol Res*, Sep;46(10):964-976 (2016) doi: 10.1111/hepr.12654.
- 6.
7. Okita K, Mizuguchi T, Ota S, Ishi M, Nishidate T, Ueki T, Meguro M, Kimura Y, Tanimizu N, Ichinohe N, Torigoe T, Kojima T, Mitaka T, Sato N, Sawada N, Hirata K. Pancreatic regeneration: basic research and gene regulation. *Surg Today*, Jun;46(6):633-640 (2016) doi: 10.1007/s00595-015-1215-2.
8. 市戸義久・石井雅之・谷水直樹・水口徹・平田公一・三高俊広. 組織幹細胞と組織修復 *Surgery Frontier*. Vol 21. No. 3, 67-69・2014
9. Tanimizu N, Mitaka T. Reevaluation of cell supply from liver stem/progenitor cells and the lineage conversion during liver development and regeneration. *Organogenesis*, 10(2): 208-215. (2014) http://dx.doi.org/10.4161/org.27591

原著論文

1. Kino J, Ichinohe N, Ishii M, Mitaka T. Isolation and expansion of rat hepatocytic progenitor cells. **Methods Mol Biol**, 1905: 29-41 (2019) doi: 10.1007/978-1-4939-8961-4_4
2. Tanimizu N. Isolation of bipotential liver progenitor cells from neonatal mouse liver. **Methods Mol Biol**, 1905: 9-17 (2019) doi: 10.1007/978-1-4939-8961-4_2.
3. Tanimizu N, Ichinohe N, Mitaka T. Intrahepatic bile ducts guide formation of the hepatic nervous system in developing mouse liver. *Development*, Apr25; 145 (9) pii: dev159095 (2018) doi: 10.1242/dev.159095
4. **Ishii M, Kino J**, Ichinohe N, Tanimizu N, Ninomiya T, Suzuki H, Mizuguchi T, Hirata K, Mitaka T. Hepatocytic parental progenitor cells of rat small hepatocytes maintain a self-renewal capability after long-term culture. *Scientific Reports*, Apr 11; 7: 46177 (2017) doi: 10.1038/srep46177
5. Ichinohe N, Ishii M, Tanimizu N, Kon J, Yoshioka Y, Ochiya T, Mizuguchi T, Hirata K, Mitaka T*. Transplantation of Thy1⁺ cells accelerates liver regeneration by enhancing the growth of small hepatocyte-like progenitor cells via IL17RB signaling. *Stem Cells*, Apr; 35(4): 920-931 (2017) doi: 10.1002/stem.2548
6. Tanimizu N, Ichinohe N, Yamamoto M, Akiyama H, Nishikawa Y, Mitaka T. Progressive induction of liver progenitor cells in chronically injured liver. *Scientific Reports*, Jan 4; 7: 39990 (2017) doi: 10.1038/srep39990
7. Tanimizu N, Ichinohe N, Ishii M, Kino J, Mizuguchi T, Hirata K, Mitaka T. Liver progenitors isolated from adult healthy mouse liver efficiently differentiate to functional hepatocytes in vitro and repopulate liver tissue. *Stem Cells*, 34(12): 2889-2901 (2016) doi: 10.1002/stem.2457
8. Tanimizu N, Kaneko K, Itoh T, Ichinohe N, Ishii M, Mizuguchi T, Hirata K, Miyajima A, Mitaka T. Intrahepatic bile ducts are developed through formation of homogeneous continuous luminal network and its dynamic rearrangement. *Hepatology*, Jul;64(1):175-188 (2016) doi: 10.1002/hep.28521
9. Kasuya J, Sudo R, Masuda G, Mitaka T, Ikeda M, Tanishita K. Reconstruction of hepatic stellate cell-incorporated liver sinusoidal structures in small hepatocyte tri-culture using microporous membranes. *J Tissue Eng Regen Med (TERMS)*, Mar; 9(3): 247-256 (2015) doi: 10.1002/term.1630
10. Tanimizu N, Kobayashi S, Ichinohe N, Mitaka T. Downregulation of miR122 by grainyhead like-2 restricts the hepatocytic differentiation potential of adult liver progenitor cells. *Development*, Dec; 141(23): 4448-4456 (2014) doi: 10.1242/dev.113654
11. Kikkawa Y, Miwa T, Tanimizu N, Kadoya Y, Ogawa T, Katagiri F, Hozumi K,

Nomizu M, Mizuguchi T, Hirata K, Mitaka T. Soluble lutheran/basal cell adhesion molecule is detectable in plasma of hepatocellular carcinoma patients and modulates cellular interaction with laminin-511 in vitro. *Exp Cell Res*, Oct 15; 328(1): 197-206 (2014) doi: 10.1016/j.yexcr.2014.07.012

12. Tanimizu N, Nishikawa Y, Ichinohe N, Akiyama H, Mitaka T. Sry HMG box protein 9-positive (Sox9+) epithelial cell adhesion molecule-negative (EpCAM-) biphenotypic cells derived from hepatocytes are involved in mouse liver regeneration. *J Biol Chem*, Mar 14; 289(11): 7589-7598 (2014) doi: 10.1074/jbc.M113.517243
13. Nakamura Y, Mizuguchi T, Tanimizu N, Ooe H, Ichinohe N, Hirata K, Mitaka T. Preoperative hepatocyte transplantation prevents cirrhotic rats from receiving the fatal damage by a liver resection. *Cell Transplantation*, 23(10): 1243-1254 (2014) doi: 10.3727/096368913X668649